

TRUCLEAR™ PRODUCT CATALOGUE



TABLE OF CONTENTS

- 2 TruClear™ Hysteroscopic Tissue Removal System**
- 3 TruClear™ 5C Hysteroscope**
- 4 TruClear™ Elite Hysteroscopes**
- 5 TruClear™ Tissue Removal Devices**
- 6 HysteroLux™ Fluid Management System**
- 7 Hysteroscopic Visualization with the TruClear™ System**
- 8 Tube Sets**

TRUCLEAR™ HYSTEROSCOPIC TISSUE REMOVAL SYSTEM



- Mechanical approach to remove intrauterine tissue.
- No scarring from thermal energy or energy discharge.
- No bubbles or electrical energy in the uterus reduces the risk of air or gas emboli and the risk of patient harm.
- Designed for visual D & C, evacuation of retained products of conception (RPOC), and fibroid and polyp removal.



Part #	Product Description	Unit of Measure
7209808	TruClear™ Control Unit	Each
7209820	TruClear™ Footswitch	Each
7209807	TruClear™ Handpiece	Each

TRUCLEAR™ 5C HYSTEROSCOPE



- Diagnostic and operative procedures with one scope.
- Lets you treat intrauterine abnormalities at the same time you diagnose them, eliminating the need for two separate scopes.
- Anatomically-designed distal end for gentle introduction into the cervical canal and uterine cavity, requiring little or no dilation.
- Enables access in high BMI patients with working channel length of 205 mm.
- Can be used in the office or the OR.


Part #	Product Description	Unit of Measure
72204752	TruClear™ Operative Hysteroscope 5C	Each
72204753	TruClear™ Sheath 5C	Each
72203506	TruClear™ Calibration Insert 2.9	Each
7205561	Hysteroscope Seal Caps	box of 5
72203618	TruClear™ Universal Hysteroscope Stopcock Replacement	Each

TRUCLEAR™ ELITE HYSTEROSCOPIES

TruClear™ Elite Hysteroscope Mini (6 mm outer diameter)



TruClear™ Elite Hysteroscope Plus (7.25 mm outer diameter)



Integrated continuous inflow and outflow with suction:

- Optimizes visibility in the uterine cavity.
- Eliminates floating tissue in the uterus for a clear operative field.
- The angled tip Enables easier insertion through the cervical canal and into the uterus.

The rod lens optical system delivers:

- High-definition (HD) images.
- New possibilities: the TruClear™ Elite hysteroscope mini has the smallest rod lens system currently available for mechanical intrauterine tissue removal.

Part #	Product Description	Unit of Measure
72204879	TruClear™ Elite hysteroscope plus (includes insert)	Each
72204878	TruClear™ Elite hysteroscope mini (includes insert)	Each
72205067	TruClear™ Elite working channel insert plus	Each
72205066	TruClear™ Elite working channel insert mini	Each
72205051	TruClear™ Elite hysteroscope seal	box of 5
90505001	TruClear™ Elite retaining cap	Each
90505000	TruClear™ Elite stopcock replacement set	Each
90505002	TruClear™ Elite Blue Ring Seals	Each

TRUCLEAR™ TISSUE REMOVAL DEVICES

TruClear™ Tissue Shaver Devices Deliver:

- Simultaneous tissue cutting and aspiration to reduce procedure time.
- Proprietary suction control feature to minimize fluid loss.
- A choice of four pathology-optimized disposable devices.

TruClear™
Soft Tissue Shaver Mini



TruClear™
Soft Tissue Shaver Plus



TruClear™
Dense Tissue Shaver Mini



TruClear™
Dense Tissue Shaver Plus



Part #	Product Description	Unit of Measure	Compatible hysteroscopes
72202536	TruClear™ Soft Tissue Shaver Mini	Each	TruClear™ Elite Hysteroscope Mini, TruClear™ Operative Hysteroscope 5.0 or 5C
7209509	TruClear™ Soft Tissue Shaver Plus	Each	TruClear™ Elite Hysteroscope Plus, TruClear™ Operative Hysteroscope 8.0
72204064	TruClear™ Dense Tissue Shaver Mini	Each	TruClear™ Elite Hysteroscope Mini, TruClear™ Operative Hysteroscope 5.0 or 5C
72203012	TruClear™ Dense Tissue Shaver Plus	Each	TruClear™ Elite Hysteroscope Plus, TruClear™ Operative Hysteroscope 8.0

HYSTEROLUX™ FLUID MANAGEMENT SYSTEM



- Delivers a high maximum flow rate of 800 mL per minute.
- Provides a rapid response to changes in intrauterine pressure.
- Has a high maximum pressure setting of 150 mm Hg.
- Enables continuous, real-time monitoring of fluid deficit when its scale is in use.
- Very low risk of fluid overload by continuously monitoring fluid volume deficit with ± 6 percent accuracy.
- Indicates critical deficit levels, pressure levels, and deficit changes with an audible alarm.
- A user-friendly touchscreen that simplifies setup.
- Integrated suction.
- Compatibility with a variety of diagnostic and operative hysteroscopes and resectoscopes to enable more procedural flexibility.

Part #	Product Description	Unit of Measure
72205000	Fluid Management System, Control Unit	Each
72205001	Fluid Management System, Scale	Each
72205015	HysteroLux™ Procedure Kit (Includes inflow tubing, outflow tubing and suction cannister accessory kit)	Each
72200024	3000 mL Hi-Flow Canisters	Box of 36
Country Kits – Power Cord and Instructions for use		
72205003	UK – EN	Each
72205004	Western Europe – EN, FR, NL, DE	Each
72205005	Mediterranean – EN, PT, ES, IT, EL, TR	Each
72205008	Nordics – EN, SV, DA, FI, NO	Each

HYSTEROSCOPIC VISUALIZATION WITH THE TRUCLEAR™ SYSTEM

One complete solution. Synchronized for clear image quality

Camera

- Capture and digitally process high-quality images without the need for upscaling. Get higher sensitivity and contrast, edge enhancement, and real color reproduction.

Clear image quality

- Achieve clear views of the uterine cavity with our single chip complementary metal-oxide-semiconductor (CMOS) technology.

Intuitive, easy-to-use interface

- Control light intensity from the camera head remote control with unique MIS-Bus system.

Robust design

- Avoid excessive heat development and mechanical shocks with the integrated electronic design.



Part #	Product Description	Unit of Measure
Flexivision 10 (FV10)		
95.3931	Camera control unit, including camera head straight with parfocal zoom, f = 14.25–28 mm; with IM capture, smoke reduction, and selective color enhancement	Each
LED Light Source		
05.0740 LED	200 LED light source with guide ADAP	Each
05.0090L	4.8 × 2300 mm hightrans light cable × 1	Each
Monitor		
MLW-2623C-DC	26 inch monitor	Each
Carts		
PT.9005.902	Cart with one monitor holder and variable height support arm	Each

TUBE SETS



For more information about all TruClear™ system products, contact your sales representative or scan the QR code.

Photo credit: iStock

Medtronic

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™ Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 20-weu-truclear-product-catalogue-4076373

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

Part #	Product Description	Unit of Measure
72205027	Hysteroscopic Vacuum Tube Set	Box of 10
		
72205028	Hysteroscopic Inflow Tube Set	Box of 10
		
7209823	Outflow Tube Set (Includes HysteroLux™ drape tube, HysteroLux™ handpiece tube and HysteroLux™ scope tube)	Box of 5
		
7209824	Suction canister accessory kit (Includes 4× jumper tubes and tissue trap)	Each
		
72205015	Hysteroscopic Procedure Kit (Includes inflow tube set, outflow tube set and suction canister accessory kit)	Box of 10
72203257	Single Cable Procedure Kit (Inclusive of one jumper tube and one tissue trap)	Each